

Single 2D barcode



SBS rack reader



AlteRead model	Single 2D barcode camera reader	SBS rack 2D barcode camera reader
Capture technology	Camera	Camera
Model	Single tube	SBS rack
Read and decode time	< 1 second	< 2 seconds
SBS format: 24, 48 & 96	n/a	Yes
Large format: 81, 100 & 192 racks	n/a	No
Custom rack	n/a	Yes
Detects missing tubes	n/a	Yes
Automatic read	n/a	Yes
Single tube read	Yes	No
Dimensions in cm	6 x 9 x 6	20.8 x 13.5 x 8
System requirements	Windows 7, 8, 10	Windows 7, 8, 10
Power	One USB cable	Universal power supply
Integration toolkits	SDK	SDK
Network connectivity	Remote control capability	Remote control capability
Environment	5°C to 35°C	5°C to 35°C
Cold rack compatible	CryoCover™	CryoCover™
Warranty	5 years	5 years