


AlteRead Single**AlteRead Wireless****AlteRead SBS rack****AlteRead Large rack**

	AlteRead Single	AlteRead Wireless	AlteRead SBS rack	AlteRead Large rack
Capture technology	Camera	Camera	Camera	Camera
Read and decode time	< 1 second	< 1 second	< 2 seconds	< 2 seconds
For use with	Single tubes	Single tubes	SBS rack	Large format & SBS racks
Barcode type	1D & 2D	1D & 2D	2D	2D
SBS format: 24, 48 & 96	n/a	n/a	Yes	Yes
Large format: 81, 100 & 192 racks	n/a	n/a	No	Yes
Custom rack	n/a	n/a	Yes	Yes
Detect missing tubes	n/a	n/a	Yes	Yes
Automatic read	Yes	Yes	Yes	Yes
Single tube read	Yes	Yes	No	No
Dimensions in cm	9.4 x 7.6 x 9.4	12.8 x 6.2 x 7	20.8 x 13.5 x 8	15.2 x 17.4 x 14.6
Power	USB cable	USB charging	Universal power supply	
System requirements	Windows / OS			
Integration toolkits	n/a	n/a	SDK	SDK
Network connectivity	Remote control capability			
Environment	5°C to 35°C			
Cold rack compatible	CryoCover™			
Warranty	5 years			