


SampleScan software

Barcode rack reader software




Included with all rack barcode readers

With minimal set-up, the intuitive project-based software is easy to use.

Labware flexibility

Compatible with all SBS and large format non-SBS racks from any manufacturer.



Automatic rack identification


Decode 1D barcodes on any rack side and 2D embedded rack barcodes.

Save projects

Easily switch between different tube rack formats.

Advanced settings

Barcode error handling, detection, reporting of missing tubes and select specific tubes for decoding.



Customise output file settings

Barcode data and output .csv file can be manually edited, and content defined by the user.

Integration and network connectivity

ActiveX and SDK for easy integration with LIMS software and other robotic systems. Network connectivity with TCP/IP.