


# SampleScan software

## Barcode rack reader software




### Included with all rack barcode readers

With minimal set-up, the intuitive project-based software is easy to use.

### Labware flexibility

Compatible with all SBS and large format non-SBS racks from any manufacturer.



### Automatic rack identification


Decode 1D barcodes on any rack side and 2D embedded rack barcodes.

### Save projects

Easily switch between different tube rack formats.

### Advanced settings

Barcode error handling, detection, reporting of missing tubes and select specific tubes for decoding.



### Customise output file settings

Barcode data and output .csv file can be manually edited, and content defined by the user.

### Integration and network connectivity

ActiveX and SDK for easy integration with LIMS software and other robotic systems. Network connectivity with TCP/IP.