

# Automated screw cap handling

## Walkaway systems for screw capping tasks

AltemisLab offers a comprehensive range of decapping and capping systems that can help with the time challenges screw cap tube users often experience.

With the AltemisLab range, automated screw cap decapping and capping is achieved for a variety of tube types and rack formats.

**Transforming lab technique.**

**XL Capping station**



**XL AutoCap**



**Advantage 96**



Tube throughput	1	8 or 12	96
Operation modes	Cap – Decap – Recap	Cap – Decap – Recap	Cap – Decap – Recap
Compatible with racked SBS tubes	No	Yes	Yes
Compatible with glass vials	Cap diameter > 10mm	No	No
Compatible with large volume polypropylene tubes	Cap diameter > 10mm	No	No
Fresh capping	Cap feeder	Cap carrier	Cap carrier
Integration friendly	Yes	Yes	Advantage 96 CELL
Dimensions - cm	43 x 35 x 56	29 x 65 x 32	35 x 31 x 48
System requirements	Windows 10, 8, 7 One USB port	Windows 10, 8, 7 One USB port	Windows 10, 8, 7 One USB or Serial port (Advantage 96 CELL only)