

Automated tube and vial labelling

Label your labware to your exact requirements

XL LabelPro™ reduces the time and effort spent on labelling your labware.

Easy to use software allows you to design your label to your exact requirements, including text, barcodes, logos and more.

Integrate with any XL tube handler for automatic labelling and walk away while it completes the repetitive work for you.

A better way



XL LabelPro



XL9 LabelPro



XL20 LabelPro



XL100 LabelPro



XL200 LabelPro

	XL LabelPro	XL9 LabelPro	XL20 LabelPro	XL100 LabelPro	XL200 LabelPro
Description	Standalone	Automated	Automated	Automated	Automated
Maximum tube/vial height - mm	125	SBS	SBS	100	125
Maximum tube/vial diameter - mm	25	8	8	22	22
Label capped or uncapped tubes	Yes	Yes	Yes	Yes	Yes
Label jacketed, label-free, flip-top tubes	Yes	Yes	Yes	Yes	Yes
Print text	Yes	Yes	Yes	Yes	Yes
Print 1D and 2D barcodes	Yes	Yes	Yes	Yes	Yes
SBS rack capacity	n/a	7	18	18	28