

Multi thawing station

Speed up sample processing & minimise data variability

Part No.	Product name
702003	Multi thawing station - Ambient
702004	Multi thawing station - Heated



Fast results

Highly uniform frozen sample thawing ensuring reproducible results. Thaw up to 1000 samples in 15 minutes.

Heated model

Combines a 600W heating module with 7 turbo fans. Suitable for samples with large molecule compounds & high melt carriers.

Easy to use

Pre-calibrated and only one switch operation.

Safe ambient thawing

Ambient air is circulated by 7 powerful fans for fast consistent thawing of all frozen sample types.

Flexible

Suitable for 96 well assay plates & racked tubes. Accessory racks are available for most non-SBS tube types.

Integration friendly

Features a remote I/O for simple integration with automation.