

Avidien microPro

Benchtop 96 channel pipettors

The world's smallest 96 channel semi-automated benchtop pipettor, available in two sizes, the microPro 20 and the microPro 300. Flexible multi-function pipetting; from simple to more complex liquid handling steps.

A better way

MicroPro 20



MicroPro 300



Volume range - μL	1-20	5-300
SBS rack formats	96/384	96/384
Compatible with deep well blocks	Yes	Yes
Tip type	Filtered and non-filtered	Filtered and non-filtered
Two position stage	Optional extra	Optional extra
Pipetting speed control	Yes	Yes
Precision	1-5 μL <3% CV 6-20 μL <2% CV	15-30 μL <2% CV 30-300 μL <1% CV
Software controls	iPad app	iPad app
iPad included	Yes	Yes
Dimensions - cm	19 x 21.6 x 32	19 x 21.6 x 32
Weight - kg	5.1	5.1
Power	Mains powered	Mains powered