

AlteCap™

Screw cap handling

Flexible systems for screw capping tasks

Automating screw cap decapping and capping is now easier and more affordable with AltemisLab's comprehensive range of decapping and capping systems. Accommodating all throughput requirements, the AlteCap™ decappers can help with the time challenges screw cap tube users often experience. All systems in the AlteCap™ range are compatible with a variety of different tube types and rack formats from any manufacturer.

Transforming lab technique

AlteCap™ Solo



AlteCap™ Switch



AlteCap™ Swift



	AlteCap™ Solo	AlteCap™ Switch	AlteCap™ Swift
Tube throughput	1	6, 8 or 12	24, 48 or 96
Operation modes	Cap – Decap – Recap	Cap – Decap – Recap	Cap – Decap – Recap
Type	Handheld	Handheld	Benchtop
Compatible with SBS racks	Yes	24, 48 or 96	24, 48 or 96
Tube manufacturer	Any	Any	Any
Fresh capping	Cap carrier	Cap carrier	Cap carrier
Stand included	Yes	Yes	N/A
Integration friendly	N/A	N/A	Swift CELL
Dimensions - cm	2.5 x 19	19 x 12.5 x 82	36.5 x 24.7 x 33.2
Power	USB charging	USB charging	Mains powered