



AlteRead Single**AlteRead Wireless****AlteRead Pendant****AlteRead SBS rack****AlteRead Large rack**

Capture technology	Camera	Camera	Camera	Camera	Camera
Read and decode time (seconds)	< 1	< 1	< 1	< 2	< 2
For use with	Single tubes	Single tubes	Single tubes	SBS rack	Large & SBS racks
Barcode type	1D & 2D	1D & 2D	1D & 2D	2D	2D
SBS format: 24, 48 & 96	n/a	n/a	n/a	Yes	Yes
Large format: 81, 100 & 196 racks	n/a	n/a	n/a	No	Yes
Custom rack	n/a	n/a	n/a	Yes	Yes
Detect missing tubes	n/a	n/a	n/a	Yes	Yes
Automatic read	Yes	Yes	Yes	Yes	Yes
Wireless receiver	No	Yes	Yes	No	No
Bluetooth connectivity	No	Yes	Yes	No	No
On board storage	No	No	Yes	No	No
Direct data export	Yes	Yes	Yes	No	No
LCD display	No	No	Yes	No	No
Dimensions in cm	9.4 x 7.6 x 9.4	12.8 x 6.2 x 7	12.9 x 5.5 x 2.7	20.8 x 13.5 x 8	15.2 x 17.4 x 14.6
Power	USB cable	USB charging	USB charging	Universal power supply	
Compatible operating systems	Windows / Mac OS		Windows / Mac OS		Windows
Integration toolkits	n/a	n/a	n/a	API	API
Network connectivity	Remote control capability				
Environment	5°C to 35°C				
Cold rack compatible	CryoCover™				
Warranty	5 years	2 years	2 years	5 years	5 years