


AlteTube™

Screw cap tubes in 96 SBS rack format



	0.5ml Screw cap	0.75ml Screw cap	0.75ml Screw cap	1.1ml Screw cap	1.4ml Screw cap	1.4ml Screw cap
Thread type	Internal	Internal	External	Internal	Internal	External
Max working volume - ml	0.3	0.49	0.55	0.75	0.92	0.95
Min storage temperature - °C	-196	-196	-196	-196	-196	-196
Tube height - mm	20.0	27	29.5	36.5	44.0	45.6
Tube height with cap - mm	27.0	34.0	32.5	43.5	51.0	48.6
Inner diameter - mm	6.8	6.8	6.5	6.8	6.8	6.5
Outer diameter with cap - mm	8.4	8.4	8.5	8.4	8.4	8.5
Tube centre to centre - mm	9.0	9.0	9.0	9.0	9.0	9.0
Identification	2D DataMatrix barcode on base					
2D DataMatrix code	White on black, 3mm x 3mm, 2 alpha prefix (AL) & 8 numeric Custom prefix available / Black on white available					
Capping options	Screw caps					
Sterility options	AlteClean™/Sterile (E-Beam)/ Ethylene Oxide (EtO) on request					
Rack format	8 x 12 - SBS (96 tubes)					
Rack material	Polypropylene base with polycarbonate lid					
Rack height with lid - mm	31	38	37	50	55.5	51.5
2D DataMatrix barcode	Bulk packed tubes, no cap – 960/case Bulk packed tubes, precapped - 960/case Racked tubes, no cap – 96/rack, 10 racks/case Racked tubes, precapped - 96/rack, 10 racks/case Bulk screw caps - 960/case Empty racks – 10 racks/case Cap carriers - 10 carriers/case					

Warning : Caps should not be immersed in Liquid Nitrogen. AlteTubes should be stored in gas phase Liquid Nitrogen (GPLN) only. For research use only.