

AlteTube™

Septum cap tubes in 96 SBS rack format

					
	0.5ml Septum cap	0.75ml Septum cap	0.75ml Septum cap [only]	1.1ml Septum cap	1.4ml Septum cap
Thread type	Internal	Internal	N/A	Internal	Internal
Max working volume - ml	0.35	0.59	0.65	0.85	0.99
Min storage temperature - °C	-80	-80	-80	-80	-80
Tube height - mm	20.0	27	33.5	36.5	44
Tube height with cap - mm	21.0	28	34.5	37.5	45
Inner diameter - mm	6.8	6.8	6.8	6.8	6.8
Outer diameter with cap - mm	8.4	8.4	8.4	8.4	8.4
Tube centre to centre - mm	9.0	9.0	9.0	9.0	9.0
Identification	2D DataMatrix barcode on base				
2D DataMatrix code	White on black, 3mm x 3mm, 2 alpha prefix (AL) & 8 numeric Custom prefix available / Black on white available				
Capping options	Septum caps				
Sterility options	AlteClean™/Sterile (E-Beam)/ Ethylene Oxide (EtO) on request				
Rack format	8 x 12 - SBS (96 tubes)				
Rack material	Polypropylene base with polycarbonate lid				
Rack height with lid - mm	26	33	39	44	50
2D DataMatrix barcode	Bulk packed tubes, no cap – 960/case Racked tubes, no cap – 96/rack, 10 racks/case Septum caps – 10 mats/case Empty racks – 10 racks/case				

For research use only.