



AlteTools™ CoboTube

Benchtop automated tube handler

Part No.	Product name
000000	AlteTools™ CoboTube

The first benchtop robot to handle the cold, capable of picking and placing tubes with precision and flexibility across a wide variety of rack formats



Enhanced productivity

Works with 4 SBS racks or Cryoracks at a time, with additional set downs available for increased capacity.

Sample Integrity

Ensure sample integrity by minimising temperature impact between samples of interest, and destination

Intuitive User Interface

Tubes to be picked can be simply selected either through the app or directly on a pick-list in Excel or csv. The robot will pick the tubes selected & place them into a destination rack in a pattern also configurable through the interface.

Integration Friendly

The System can be integrated into a laboratory workflow primarily for handling sample tubes & is equipped to work seamlessly with cryogenic or dry ice containers.

Tube Sorting

Enables end users to reformat, cherry-pick or re-array directly from SBS 24, 48 & 96-well racks, & 9x9 or 10x10 Cryoracks, with ease.

Modular Product

With options to incorporate single tube readers, rack readers, or integrate a weigh station, the unit can be customised to align with your workflow.

Flexible

Easily customisable to optimise the platform size for any workflow.