

AlteRead™ Single 2D barcode readers

AlteRead™ Single



AlteRead™ Wireless



AlteRead™ Pendant



Capture technology	Camera	Camera	Camera
Read and decode time	< 1 second	< 1 second	< 1 seconds
For use with	Barcodes on tubes & racks	Barcodes on tubes & racks	Barcodes on tubes & racks
Barcode type	1D & 2D	1D & 2D	1D & 2D
Wireless Receiver	n/a	Yes	Yes
Bluetooth connectivity	n/a	Yes	Yes
IR sensing technology	Yes	Yes	Yes
Direct data export	Yes	Yes	Yes
On board storage	No	No	Yes
LCD display	No	No	Yes
Dimensions in cm	9.4 x 7.6 x 9.4	12.8 x 6.2 x 7	12.9 x 5.5 x 2.7
Power	USB cable	USB charging	USB charging
Compatible operating systems	Windows / Mac OS	Windows / Mac OS / Google Android / Apple iOS	
Warranty	5 years	2 years	
Environment	5°C to 35°C		
Cold rack compatible	CryoCover™		

AlteRead™ Pendant

Single tube 2D & 1D barcode camera reader

Part No.	Product name
255002	AlteRead™ Pendant Single tube 2D & 1D barcode camera reader



AlteRead™ pendant

The first totally portable single 2D & 1D barcode reader.

Pocket size

Small and portable, with convenient lanyard included.

LCD Display

Shows 1D & 2D barcodes instantly on the easy to read screen.

Check barcodes anywhere

Read barcodes anywhere, no need for a computer connection.

IR sensing technology

Tubes and racks are instantly recognised for decoding.

Fast decoding

Decode a single tube 2D or 1D barcode in less than 1 second.

On-board storage

Option to store 2000 barcodes and download the results at a later time.

Plug & Play


No software to install and display results in any application.

Versatile

Connect to computers, phones and smart devices via Bluetooth.

Tube compatibility

All 2D and 1D barcoded tubes from any manufacturer.



See the AlteRead™ range of 2D & 1D barcode readers in the AltemisLab product catalogue